


APPLICATION SHEET

Glass / Ceramics – DMA EPLEXOR®

Investigations on Fused Silica

With the EPLEXOR® HT series, NETZSCH GABO Instruments GmbH shows worldwide the first DMA instrument providing a temperature range up to 1500°C. This class of instruments is aimed for providing investigations on glass, metals, ceramics and even high-temperature composite materials in the high-temperature range. The figure shows the test


results obtained for different fused silica. Rectangular glass samples with dimensions of 5 mm (width) x 2 mm (thickness) x 30 mm (length) were investigated. For both samples, a test frequency of 1 Hz and a heating rate of 2 K/min were applied. The sample holder used for this test was a newly designed asymmetrical 3-point bending device. The damping properties ($\tan\delta$) as a function of temperature significantly differ from each other.



Resolution – With Which Accuracy Can $\tan\delta$ Be Determined?

This picture shows storage modulus E' and loss factor $\tan\delta$ of a glass, measured in the 3-point bending mode. The experiment provides the low absolute values for the visco-elastic damping as a function of temperature.

Such minimal $\tan\delta$ values, corresponding to a phase shift between force and sample deformation, can only be determined with an extremely high resolution in the phase shift measurement. We therefore have equipped our EPLEXOR® series with high precision phase angle determination. The resolution is better than $\tan\delta \leq 10^{-4}$.



APPLICATION SHEET

Glass / Ceramics – DMA EPLEXOR®

A New Tool for Tests on Very Stiff, but Brittle Low-Damping Materials – Asymmetric 3-Point Bending Device

The picture shows the complex modulus E^* and loss factor $\tan\delta$ of a glass measured in the asymmetric 3-point bending mode. The observed $\tan\delta$ peak with a maximum at 200°C (below the glass transition of over 400°C) has to be

related to thermally activated molecular motions of the Na^+ ions within the glass.

Friction processes due to the interaction between the sample and sample holder can totally be suppressed. Therefore, the asymmetric 3-point bending device enables the user to investigate brittle low-damping materials ($\tan\delta \sim 0.0002$) with very high accuracy.

