

**NETZSCH**

Proven Excellence.



## TMA 402 **F3** Hyperion<sup>®</sup> Polymer Edition

Thermomechanical Analysis – TMA  
Tailor-Made for Low-Temperature Applications

Analyzing & Testing

# Tailored to Polymer Applications

For research and development,  
simulation input and quality control.

## TMA Tests on Polymers

- Coefficient of linear thermal expansion
- Phase transition temperatures
- Glass transition temperatures
- Dilatometric softening point
- Volumetric expansion
- Density changes
- Delamination
- Isostrain
- Creep
- Relaxation
- Stress/strain curve
- Melting and crystallization

Thermomechanical analysis (TMA) is a technique for determining the dimensional changes in solids, liquids or pasty materials as a function of temperature and/or time under a defined mechanical force (DIN 51005, ASTM E 831, ASTM D696, ASTM D3386, ISO 11359 – Parts 1 to 3). TMA has proven to be the preferred method for measuring the coefficient of linear thermal expansion in plastics, a fundamental material property describing the ability of a polymer to expand when exposed to varying temperatures.

## TMA Analysis for Polymers


Polymers undergo changes in their thermomechanical properties during heating and cooling. TMA analyses can give insights into molecular orientation and quenching effects during cooling, thus allowing the design of adhesives and other hybrid joints as well as for quality control of shrink films. Such analyses can provide valuable insight into the composition, structure, production conditions and application possibilities for various materials.

When it comes to the detection of  $T_g$  for polymers with fillers or highly crosslinked materials, such as composites or printed circuit boards, TMA offers a higher degree of sensitivity than the DSC method. TMA is therefore a perfect method to be used to detect weak physical transitions that are associated with changes in modulus, post-curing and delamination, which are hard or impossible to measure with other methods.

## Thermal Expansion

The linear thermal expansion is an important variable for assessing the dimensional change of a material in response to a change in temperature. It shows how much a material will shrink or expand during processing, whether dissimilar materials can be joined, where the phase change occurs and where the CTE changes.


This plot shows the thermal expansion ( $dL/L_0$  in %) of an epoxy resin between  $-70^\circ\text{C}$  and  $270^\circ\text{C}$ . In the first heating (blue curve), the onset of the  $T_g$  occurs at  $123^\circ\text{C}$ . In the second heating (red curve), the onset of  $T_g$  is slightly shifted, to  $125^\circ\text{C}$ . This shift could be due to relaxation effects or post-curing.



CTE measurement on an epoxy resin with a sample length of 6 mm in expansion mode (fused silica sample holder); 1<sup>st</sup> and 2<sup>nd</sup> heating runs at a rate of 2 K/min

## TMA – THE METHOD PRECISELY DETERMINES DIMENSIONAL CHANGES OF PLASTICS, ELASTOMERS AND COMPOSITE MATERIALS

### OPERATING PRINCIPLE



Irrespective of the type of deformation selected (expansion, compression, penetration, tension or bending), every length change in the sample is communicated to a highly sensitive inductive displacement transducer (LVDT) via a pushrod and transformed into a digital signal.

The pushrod and corresponding fused silica sample holders can be quickly and easily changed out, in order to optimize the system for the respective application.

# TMA 402 *F3 Hyperion*<sup>®</sup> Polymer Edition

## Gaining Valuable Information about

### Detect Even the Slightest Dimensional Changes

The LVDT constitutes the centerpiece of the NETZSCH TMA 402 *F3 Hyperion*<sup>®</sup> Polymer Edition. The technology behind it is tried-and-true: Even the slightest of length changes, into the nanometer range (digital resolution of 0.125 nm), can be measured and detected.

### No-Hassle Cooling to -70°C

The TMA 402 *F3 Hyperion*<sup>®</sup> Polymer Edition is specifically designed for polymer applications. It comes with a compact, highly reactive furnace capable of covering a temperature range from -70°C to 450°C and uses a mechanical cooling system without the need for LN<sub>2</sub>.

### Flexible Atmospheres in a Vacuum-Tight TMA System

All joints are designed to be vacuum-tight, allowing for measurements in a highly pure atmosphere or under vacuum. This TMA includes one software-controlled mass flow controller (MFC) and can be upgraded with another independent MFC. This provides optimum flexibility in gas control and in changing the purge and atmospheric gas rate.

Made in Germany  
by NETZSCH




# Product Performance and Processing Behavior

## Determination of Visco-Elastic Properties like Relaxation, Creep and Stress/Strain

The TMA 402 **F3** Hyperion® Polymer Edition now offers not only the ability to keep the force constant and to measure the length change, but also to change the dL displacement and measure the corresponding force. This can, for example, be used in a stress relaxation test where a sample is stretched by a specific amount at a defined temperature. During the test, the deformation is kept constant and the progression of the force is recorded. This force continuously decreases as a result of material relaxation. The stress-relaxation strength is ultimately defined by the residual stress measured after a defined exposure period. The data can be depicted graphically in a stress-time diagram. It is then possible to read off both the stress-relaxation behavior and the values for the relaxation rate and time.

## Precise Force Control Enables Tests on Sensitive Materials


The force operating on the sample is generated electromagnetically. This ensures a quick response time for experiments with a changing load. A highly sensitive force sensor (digital resolution < 0.01 mN, max. force ±3 N) continuously measures the force exerted via the pushrod and readjusts it automatically. This sets the TMA 402 **F3** Hyperion® Polymer Edition apart from other instruments which only use preset values. The control system allows users to set the force value in the mN-range. This enables testing even on sensitive materials such as thin fibers or films. The force being exerted upon the sample can be altered via the software in a stepwise or linear fashion. This makes it especially simple to carry out tests such as creep.



# Accessories for All Application Needs

## Easy-to-Change Sample Holders for Various Sample Geometries and Measuring Modes

Fused silica sample holders suitable for the TMA 402 F3 Hyperion® Polymer Edition



Expansion/  
Compression  
– pushrod  
with flat tip,  
Ø 4 mm

Penetration –  
pushrod with flat tip,  
Ø 1 mm and Ø 4 mm  
hemispherical tip

3-Point  
Bending,  
5 mm

Tension

The **expansion/compression** mode is used for samples with different geometries, such as cylindrical or rectangular.

The **penetration** mode determines the softening point of a sample. For this test, a hemispherical or cylindrical probe can be used.

**3-point bending** sample holder with a bending length of 5 mm is available.

The **tension** mode is used to measure expansion and shrinkage, creep and relaxation behavior on thin films or fibers.

**New tension sample holder with fixed lower clamp for almost force-free tension measurements**



Sample alignment tool for thin and sensitive samples



New tension sample holder




# Defining Ideal Production and Service Conditions for Polymer Materials

## APPLICATIONS

### Relaxation – An Important Material Property for Polymer Films

Plastic packaging needs to be flexible, lightweight, strong, sometimes impermeable and easy to sterilize. Both the material used and the process conditions determine the resulting property portfolio. It is important to have knowledge about a film material's relaxation properties; this allows for predictions regarding its fatigue and wear resistance during use. Shrink films apply a film's tendency to relax at elevated temperatures after having been pre-stretched during processing. Over time, shrink wrap tends to loosen due to creep and stress relaxation. In the test shown here, the sample was kept under constant strain and the progression of the tensile force was recorded.



Relaxation measurement on LDPE film by use of displacement step, fused silica tension sample holder. Measured at room temperature under a nitrogen atmosphere; sample length 10 mm; sample thickness 5 mm.






## Anisotropic Behavior in Composite Systems

TMA can be used to characterize multiple properties of composite systems. It can be used to determine the glass transition temperature ( $T_g$ ) in a thermoset or thermoplastic matrix composite. Determining the start of the glass transition is a reliable method for determining the upper working temperature limit of a thermoset material. For thermoplastic matrix composites such as PP-GF,  $T_g$  signifies the temperature at which the material starts to behave in a more ductile fashion.

Fibers and other fillers significantly reduce thermal expansion. The degree of anisotropy of the filler and the filler orientation both have a great impact on the linear coefficient of thermal expansion. For example, uni-directional composites exhibit the CTE of the fibers in the fiber direction. They also exhibit a mixture of the matrix and fiber CTE as a function of fiber




Measurement on a PP-GF-UD composite material. Sample size 25 mm, heating rate 5 K/min from -70°C to 140°C, N<sub>2</sub> atmosphere, expansion sample holder made of fused silica.

volume content perpendicular to the fiber direction. The red curve depicts the measurement in the fiber direction; the CTE is dominated by the low thermal expansion of glass. The sample

that is 90° to the fiber direction is dominated by the polypropylene matrix and shows a much higher CTE. Therefore, it is only possible to observe the  $T_g$  of PP in this measurement.



## Specifying the Right Material for Functional Electronics



Determination of time to delamination on an FR4 composite circuit board. Sample size 6.35 mm<sup>2</sup> as defined by IPC standard, dried for 2 hours pre-measurement at 105°C, heating rate 10 K/min, N<sub>2</sub> atmosphere, sample holder made of fused silica. Ramped to isotherm at 260°C and 300°C. The green curve at 260°C (IPC standard) shows no delamination. At 300°C, delamination starts after 18.1 min (red curve).

**Thermal expansion is a major factor in the failure of electronic products.**

Industry standards require the measurement of thermal expansion, glass transition and softening point under IPC standards [see IPC-TM-650 2.4.24.1 Time to Delamination (TMA Method)]. After the changeover to lead-free soldering processes in circuit board production, the melting temperature of soldering materials increased. This led to the delamination of printed circuit boards and assemblies due to the higher thermal load. Manufacturers had to react by changing to FR4 substrates with a higher  $T_g$ . Nevertheless, even today most FR4 substrates are ordered with general material specifications and can vary in their material


properties. At the glass transition event, the rate of expansion for the epoxy matrix increases which can lead to delamination between the fibers and matrix and consequently to product failure.

The figure shows a measurement on an FR4 composite in which the time to delamination was recorded. Two measurements were conducted. One with an isothermal temperature of 260°C and a second one with an isothermal temperature of 300°C. In the first measurement, no delamination effect was visible. In the second measurement, the time to delamination was recorded at 18.1 min after being held at an isothermal temperature of 300°C.



Typical de-coloration as a sign of beginning degradation: left: before measurement, middle: after being measured at 260°C, right: after being measured at 300°C. Neither of the measured samples shows any visible delamination, although the TMA method is sensitive enough to detect it at 300°C.

## Test on Heat-Shrink Tubing



Displacement control on heat-shrink tubing using a tension sample holder made of fused silica, with displacement kept constant while measuring force. Samples were heated from RT to 120°C, atmosphere N<sub>2</sub>, sample size 25 mm (blue) and 10 mm (green).


Heat-shrink tubing, also known as shrink sleeve, is used to repair and insulate wires and cables. After sliding the tubing onto the cable, a heat source is used to make it shrink and create a tight seal. Heat-shrink tubing, by its very nature, is stretchy and changes its shape. A TMA can help to gather information about the temperature at which the material starts to shrink, how much it shrinks and with how much force. The plot shows the two samples measured. The first sample (green) starts to shrink at around 60°C with a shrinkage of 40% by the end of the measurement. In the second sample, the displacement was kept constant and the corresponding force measured. A maximum force of 0.976 N was recorded.





## Thermal Expansion of Insulation Material

Polymer foams are used where light weight is required along with an excellent strength-to-weight ratio, superior thermal and acoustical insulation properties, and good energy absorption behavior. Expanded foams such as EPS and EPP are gaining in importance not only in the building, but also in the automotive industry. The thermal expansion is an important quality criterion when choosing insulation material, since expansion behavior differs significantly among materials with different densities. The foam with the lowest density shows the highest CTE while still having a similar  $T_g$ .



Expansion measurement on three different EPP foams with different densities. Sample size 20 mm, heating rate: 5 K/min starting from -40°C to 80°C, mechanical cooling under  $N_2$  atmosphere, expansion sample holder made of fused silica and 50 mN sample load.



# Proteus® – Only One Click Away From

## Input Assistant for a Fast Measurement Start and Method-Based Automatic Evaluation

The *Proteus*® software allows for properties and methods from previously executed measurement files to be applied with a simple mouse click. The evaluation steps for a reference test run can be saved in a method and applied, fully automatically, to a sample measurement after its termination. It is also possible to have the software highlight any results deviating from the selected quality criteria.

## AutoEvaluation – The One-Click Evaluation

*AutoEvaluation* is an intelligent software functionality exclusively offered by NETZSCH. It is a self-acting evaluation of thermoanalytical measurement curves that works without using pre-defined macros. This is an immense support and time saver. *AutoEvaluation* offers special functions for the evaluation of various materials. When measuring polymers, *AutoEvaluation* will automatically find the onset of the glass transition as well as the peaks of softening points and evaluate these with just one click.

## Temperature-Modulated TMA

For temperature-modulated TMA measurements, the modulation amplitude and period can be defined segment by segment. The evaluation software allows for the determination of total TMA, reversing TMA, non-reversing TMA, total CTE, reversing CTE, non-reversing CTE, amplitude and phase with graphic display of the results curve in multi-window technique. *Proteus*® also offers the ability to export graphs and print out or export data as ASCII files.

### At a Glance – Highlights of the TMA *Proteus*® Software TMA 402 **F3** Hyperion® Polymer Edition

Automatic sample length detection	■
Force adjustment/segment	■
Softening point detection	■
<i>c-DTA</i> ®	□
Force modulation	□
Temperature modulation	□
Strain control	■
Report generator	■
<i>Identify</i>	□
<i>AutoEvaluation</i>	■

■ Included in standard configuration  
□ Optional

# CLEVER FEATURES FOR INTELLIGENT ANALYSIS

## Reliable Results

### Identify – Identification and Classification of TMA Curves

The *Identify* database offers a state-of-the-art means of verifying materials; it allows for the comparison of a given curve to other individual curves (e.g., groups of curves in quality control) or to literature data from selected libraries. Any libraries and classes created by the user can be edited or expanded within *Identify*.

Database entries can be filtered by a variety of criteria, and measurement curves – even those of different types – can be overlapped for purposes of comparison.

### Software Options for Advanced Evaluation Steps

*Kinetics Neo* is a software option offered by NETZSCH for the characterization of curing and degradation behavior.

Automatic determination of initial sample length in expansion, penetration and tension modes!

### TMA 402 F3 Hyperion® Polymer Edition

Furnace	IC-furnace: -70°C to 450°C using mechanical cooling
Heating rates	0.001 K/min to 30 K/min
Measuring ranges/ $\Delta l$ resolution	▪ 500 $\mu\text{m}$ / 0.125 nm ▪ 5000 $\mu\text{m}$ ( $\pm 2500 \mu\text{m}$ ) / 1.25 nm
Force and displacement	Simultaneous measurement of force and displacement signal
Force range (load at sample)	0.001 N to 3 N without using additional weights
Force resolution	< 0.01 mN
Interchangeable sample holder systems	Fused silica, up to 450°C
Sample dimensions	▪ Expansion/penetration: Length: 30 mm max.; sample holder $\varnothing$ 8 mm ▪ Tension: Length: 30 mm max.; width: 6 mm, thickness: 1 mm ▪ 3-point bending: Length: 10 mm max.; width: 5 mm
Atmospheres	1 mass flow controller (MFC, 2 <sup>nd</sup> independent MFC optional), gas flow rate 0 to 250 ml/min (software-controlled), inert, oxidizing, reducing, vacuum ( $10^{-4}$ mbar)
Various accessories	Spacers, crucibles and special wax containers





# Expertise in Service

## Our Expertise – Service

All over the world, the name NETZSCH stands for comprehensive support and reliable service, before and after sale. Our qualified personnel from the technical service and application departments are always available for consultation.

In special training programs tailored for you and your employees, you will learn to tap the full potential of your instrument.

To maintain and protect your investment, you will be accompanied by our experienced service team over the entire life span of your instrument.

## Our Expertise – Applications Laboratories

The NETZSCH Thermal Analysis applications laboratories are a proficient partner for nearly any Thermal Analysis issue. Our involvement in your projects begins with proper sample preparation and continues through meticulous examination and interpretation of the measurement results. Our diverse methods and over 30 different state-of-the-art measuring stations will provide ready-made solutions for all your thermal needs.

Within the realm of thermal analysis and the measurement of thermophysical properties, we offer you a comprehensive line of the most diverse analysis techniques for materials characterization.

Measurements can be carried out on samples of the most varied of geometries and configurations. You will receive high-precision measurement results and valuable interpretations from us in the shortest possible time. This will enable you to precisely characterize new materials and components before actual deployment, minimize risks of failure, and gain decisive advantages over your competitors.

### TECHNICAL SERVICE



Maintenance and Repair



Software Updates



Exchange Service



IQ/OQ Documents



Calibration Service



Spare Part Assistance



Moving Service

### TRAINING



Training




Comprehensive Instrument and Method Training

### LABORATORY



Application Service and Contract Testing



The NETZSCH Group is an owner-managed, international technology company with headquarters in Germany. The Business Units Analyzing & Testing, Grinding & Dispersing and Pumps & Systems represent customized solutions at the highest level. More than 3,700 employees in 36 countries and a worldwide sales and service network ensure customer proximity and competent service.

Our performance standards are high. We promise our customers Proven Excellence – exceptional performance in everything we do, proven time and again since 1873.

When it comes to Thermal Analysis, Calorimetry (adiabatic & reaction), the determination of Thermophysical Properties, Rheology and Fire Testing, NETZSCH has it covered. Our 50 years of applications experience, broad state-of-the-art product line and comprehensive service offerings ensure that our solutions will not only meet your every requirement but also exceed your every expectation.

## Proven Excellence.■

NETZSCH-Gerätebau GmbH  
Wittelsbacherstraße 42  
95100 Selb  
Germany  
Tel.: +49 9287 881-0  
Fax: +49 9287 881 505  
at@netsch.com

**NETZSCH**®

[www.netsch.com](http://www.netsch.com)