

Technical Specifications

Plate/Sheet Size Supported	Optimized for two plate formats 20 x 20 cm and 20 x 10 cm
Plate/Sheet Material Supported	Glass, aluminium and plastic
Derivatization method	Automatic spraying using unique “micro droplet” technology
Reagent Consumption	4 ml for 20x20cm plates 2-3 ml for 20x10cm plates
Instrument operation and control	Via nine illuminated buttons (stand-alone instrument)
Operating temperature	10–40°C
Power connection	100–240 VAC, 50–60 Hz, 30 W
Dimensions (W x D x H)	245 x 430 x 355 mm
Weight	12 kg